

Product Name	Recombinant Human CD40	SDS-PAGE gel	
Synonym(s)	Tumor necrosis factor receptor superfamily member 5; TNFRSF5	1	2
Quantity	50 µg, 100 µg	kDa	M
Species	Human	116	
Class Type	Recombinant Protein	66.2	
Molecular weight	~32 kDa (due to glycosylation), 20.7 kDa (predicted)	45.0	
Purity	>95% by SDS-PAGE and HPLC	35.0	
Tag	C-terminal Polyhistidine Tag	25.0	
Expression Source	HEK 293 cells	18.4	
GenBank Accession #	NP_001241.1, a. a. 1-194	14.4	
Application	Western Blotting, Functional Studies. ELISA		
Formulation	Lyophilized from sterile PBS, pH 7.4		
Storage and Stability	Stable for 12 months at -80°C, Avoid freeze/thaw cycles	1 – MW Marker	2 – CD40
Description	CD40 is a receptor protein found on the surface of immune cells, including B cells, dendritic cells, and macrophages. It plays an important role in regulating immune responses, particularly in the activation and proliferation of B cells. When CD40 is activated by its ligand, CD154, it leads to the upregulation of costimulatory molecules and cytokine production, promoting an effective immune response. Dysregulation of the CD40/CD154 pathway has been implicated in various autoimmune and inflammatory disorders, as well as cancer. Targeting this pathway has therapeutic potential for these diseases.		
Reference	Bretscher, P A. <i>Proceedings of the National Academy of Sciences of the United States of America</i> vol. 96,1 (1999): 185-90. Wennhold, Kerstin et al. <i>Oncotarget</i> vol. 8,17 (2017): 27740-27753. Elgueta, Raul et al. <i>Immunological reviews</i> vol. 229,1 (2009): 152-72. Vonderheide, Robert H. <i>Annual review of medicine</i> vol. 71 (2020): 47-58.		

Related products:

Product Name	Catalog #	Size
Recombinant Human PD-L1	237351	100 ug
Recombinant Human CD40	232340	100 ug
Recombinant Human CD40L	2323405	100 ug
Recombinant Human CD70	232370	100 ug
Recombinant Human CD155	2323155	100 ug
Recombinant Human GITRL	2344875	100 ug
Recombinant Human OX40	236940	100 ug
Recombinant Human OX40L	2369405	100 ug

DNA Polymerase Theta Activity Assay Kit	362101	96 reactions
DNA Polymerase Theta-N-Helicase Domain	7657643	20 ug, 100 ug
DNA Polymerase Theta-C terminal Domain	7657283	20 ug, 100 ug, 1 mg
DNA Polymerase Theta Full Length protein	7657385	10 ug, 50 ug
T7 RNA polymerase	777627	5000 U, 25000U, 100000U
T7 High Yield RNA Synthesis Kit	777627-RK	25 rxns, 50 rxns, 100 rxns
Recombinant Mouse Leukemia Inhibitory Factor (mLIF)	11-0001	10 ug, 100 ug
Recombinant Human LIF	12-0002	10 ug, 100 ug, 1 mg
Recombinant Human LIF, Animal-Free	12-0002AFR	10 ug, 100 ug, 1 mg
Recombinant Human FGF-basic, Carrier-free	12-0005CFR	50 ug, 100 ug, 500 ug, 1 mg
Kras Wild Type (WT), GST-tag	5727-4121G	50 µg, 100 µg
Kras WT, GST-tag, GDP Loaded	5727-WTG-G	50 µg, 100 µg
Kras WT, GST-tag, GppNHp loaded	5727-WTG-GP	50 µg, 100 µg
Kras G12C, GST-tag	5727-4122G	50 µg, 100 µg
Kras G12C, GST-tag, GDP Loaded	5727-4122G -G	50 µg, 100 µg
Kras G12C, GST-tag, GppNHp loaded	5727-4122G -GP	50 µg, 100 µg
Kras G12C Nucleotide Exchange Assay Kit	5727-4122NK	384 reactions
Kras G12C – cRAF Binding Assay Kit	5727-4122BK	384 reactions
Human RBD-RAF1, N-His tag, C-FLAG tag	7237231	50 µg, 100 µg
Human SOS1, No Tag	7671	50 µg, 100 µg
Human SOS1, His-Avi-Tag	7671HA	50 µg, 100 µg
Recombinant SARS-CoV-2 Mpro, 3CL protease	728201	50 ug, 500 ug
SARS-CoV-2 Mpro (3CL Protease) Assay Kit	728203	96 reactions

This product is for research use only and not for diagnostic or therapeutic use.