



Product Name	Recombinant Human Kras Wild-Type (WT), GST-Tag, GppNHp loaded	<p>SDS-PAGE gel</p>  <p>M – MW Marker 1 – GST-Kras WT</p>
Synonym(s)	KRAS2, RASK2, C-K-RAS, CFC2, KRAS proto-oncogene	
Quantity	50 µg, 100 µg	
GenBank #	NM_004985.5	
UniProt #	P01116	
Amino Acids	a. a. 1-169	
Molecular weight	45.94 kDa + GppNHp	
Purity	>95% by SDS-PAGE	
Tag	N-GST Tagged	
Expression Source	E. coli	
Application	Inhibitor Screening, binding assays, enzyme kinetics, and selectivity profiling.	
Formulation	20 mM HEPES, 50 mM NaCl, 1 mM DTT, 10% glycerol, pH 7.5	
Storage and Stability	Stable for at least 6 months at -80°C, Avoid freeze/thaw cycles	
Description	Kras is a member of the RAS protein family, which are a class of small GTPases involved in cell signaling. Plays a role in promoting oncogenic events by inducing transcriptional silencing of tumor suppressor genes (TSGs) in colorectal cancer (CRC) cells in a ZNF304-dependent manner	
Reference	<p>Roman C. Hillig, et al. PNAS. 2019, vol. 116, no. 7 , 2551–2560 Kari Kopra, et al. Anal. Chem. 2020, 92, 4971–4979 Fernández-Medarde A., et al. Genes (Basel). 2021;12(5):681.</p>	

Kras (WT)-cRAF Binding



The TR-FRET binding assay of Kras (WT) and cRAF was performed at 10 nM cRAF and indicated amount of GDP-loaded Kras and GppNHp-loaded Kras.

Related products:

<u>Product Name</u>	<u>Catalog #</u>	<u>Size</u>
Kras Wild Type (WT), GST-tag	5727-4121G	50 µg, 100 µg
Kras WT, GST-tag, GDP Loaded	5727-WTG-G	50 µg, 100 µg
Kras WT, GST-tag, GppNHp loaded	5727-WTG-GP	50 µg, 100 µg
Kras G12C, His -tag	5727-4122H	50 µg, 100 µg
Kras G12C, GST-tag	5727-4122G	50 µg, 100 µg
Kras G12C, GST-tag, GDP Loaded	5727-4122G -G	50 µg, 100 µg
Kras G12C, GST-tag, GppNHp loaded	5727-4122G -GP	50 µg, 100 µg
Kras G12D, GST-tag	5727-4123G	50 µg, 100 µg
Kras G12D, GST-tag, GDP Loaded	5727-4123G -G	50 µg, 100 µg
Kras G12D, GST-tag, GppNHp loaded	5727-4123G -GP	50 µg, 100 µg
Kras G12R, GST-tag,	5727-4127G	50 µg, 100 µg
Kras G12R, GST-tag, GDP Loaded	5727-4127G -G	50 µg, 100 µg
Kras G12R, GST-tag, GppNHp loaded	5727-4127G -GP	50 µg, 100 µg
Kras G12V, GST-tag,	5727-4128G	50 µg, 100 µg
Kras G12V, GST-tag, GDP Loaded	5727-4128G -G	50 µg, 100 µg
Kras G12V, GST-tag, GppNHp loaded	5727-4128G -GP	50 µg, 100 µg
Kras WT Nucleotide Exchange Assay Kit	5727-4121NK	384 reactions
Kras G12C Nucleotide Exchange Assay Kit	5727-4122NK	384 reactions
Kras G12D Nucleotide Exchange Assay Kit	5727-4123NK	384 reactions
Kras G12R Nucleotide Exchange Assay Kit	5727-4127NK	384 reactions
Kras G12V Nucleotide Exchange Assay Kit	5727-4128NK	384 reactions
Kras WT – cRAF Binding Assay Kit	5727-4121BK	384 reactions
Kras G12C – cRAF Binding Assay Kit	5727-4122BK	384 reactions
Kras G12D– cRAF Binding Assay Kit	5727-4123BK	384 reactions
Kras G12R– cRAF Binding Assay Kit	5727-4127BK	384 reactions
Kras G12V – cRAF Binding Assay Kit	5727-4128BK	384 reactions
Human RBD-RAF1, N-His tag, C-FLAG tag	7237231	50 µg, 100 µg
Human SOS1, No tag	7671	50 µg, 100 µg
Human SOS1, Avi-His tag	7671HA	50 µg, 100 µg

Products are for research use only and are not intended for human use. We do not sell to patients.