

<b>Product Name</b>	Recombinant Human Kras Wild-Type (WT), GST-Tag	
<b>Synonym(s)</b>	KRAS2, RASK2, C-K-RAS, CFC2, KRAS proto-oncogene	
<b>Quantity</b>	50 µg, 100 µg	
<b>GenBank #</b>	NM_004985.5	
<b>UniProt #</b>	P01116	
<b>Amino Acids</b>	a. a. 1-169	
<b>Molecular weight</b>	45.94 kDa	
<b>Purity</b>	>95% by SDS-PAGE	
<b>Tag</b>	N-GST Tagged	
<b>Expression Source</b>	E. coli	
<b>Application</b>	Inhibitor Screening, binding assays, GDP/GTP nucleotide exchange reaction, ELISA, Immunoaffinity Purification, Protein Array, Western Blot.	
<b>Formulation</b>	PBS + 20% glycerol	
<b>Storage and Stability</b>	Stable for at least 6 months at -80°C, Avoid freeze/thaw cycles	
<b>Description</b>	Kras is a member of the RAS protein family, which are a class of small GTPases involved in cell signaling. Plays a role in promoting oncogenic events by inducing transcriptional silencing of tumor suppressor genes (TSGs) in colorectal cancer (CRC) cells in a ZNF304-dependent manner	
<b>Reference</b>	Roman C. Hillig, et al. PNAS. 2019, vol. 116, no. 7, 2551–2560 Kari Kopra, et al. Anal. Chem. 2020, 92, 4971–4979 Fernández-Medarde A., et al. Genes (Basel). 2021;12(5):681.	

**Related products:**

<b>Product Name</b>	<b>Catalog #</b>	<b>Size</b>
Kras Wild Type (WT), GST-tag	5727-4121G	50 µg, 100 µg
Kras WT, GST-tag, GDP Loaded	5727-WTG-G	50 µg, 100 µg
Kras WT, GST-tag, GppNHp loaded	5727-WTG-GP	50 µg, 100 µg
Kras G12C, His-tag	5727-4122H	50 µg, 100 µg
Kras G12C, GST-tag	5727-4122G	50 µg, 100 µg
Kras G12C, GST-tag, GDP Loaded	5727-4122G -G	50 µg, 100 µg
Kras G12C, GST-tag, GppNHp loaded	5727-4122G -GP	50 µg, 100 µg
Kras G12D, GST-tag	5727-4123G	50 µg, 100 µg
Kras G12D, GST-tag, GDP Loaded	5727-4123G -G	50 µg, 100 µg

Kras G12D, GST-tag, GppNHp loaded	5727-4123G -GP	50 µg, 100 µg
Kras G12R, GST-tag,	5727-4127G	50 µg, 100 µg
Kras G12R, GST-tag, GDP Loaded	5727-4127G -G	50 µg, 100 µg
Kras G12R, GST-tag, GppNHp loaded	5727-4127G -GP	50 µg, 100 µg
Kras G12V, GST-tag,	5727-4128G	50 µg, 100 µg
Kras G12V, GST-tag, GDP Loaded	5727-4128G -G	50 µg, 100 µg
Kras G12V, GST-tag, GppNHp loaded	5727-4128G -GP	50 µg, 100 µg
Kras WT Nucleotide Exchange Assay Kit	5727-4121NK	384 reactions
Kras G12C Nucleotide Exchange Assay Kit	5727-4122NK	384 reactions
Kras G12D Nucleotide Exchange Assay Kit	5727-4123NK	384 reactions
Kras G12V Nucleotide Exchange Assay Kit	5727-4128NK	384 reactions
Kras WT – cRAF Binding Assay Kit	5727-4121BK	384 reactions
Kras G12C – cRAF Binding Assay Kit	5727-4122BK	384 reactions
Kras G12D– cRAF Binding Assay Kit	5727-4123BK	384 reactions
Kras G12V – cRAF Binding Assay Kit	5727-4128BK	384 reactions
Human RBD-RAF1, N-His tag, C-FLAG tag	7237231	50 µg, 100 µg
Human SOS1, No Tag	7671	50 µg, 100 µg
Human SOS1, His-Avi-Tag	7671HA	50 µg, 100 µg
TEV-His	190001-R	1,000 Units, 10,000 Units
Recombinant SUMO Protease (Ulp1)	190003	1,000 units, 10,000 units
PreScission Protease (HRV 3C)	190002-R	1,000 units, 10,000 units
Recombinant YopH	200100	10 ug, 20 ug, 100 ug, 1 mg
Recombinant Mouse Leukemia Inhibitory Factor (mLIF)	11-0001	10 ug, 100 ug
Recombinant Human LIF	12-0002	10 ug, 100 ug, 1 mg
Recombinant Human FGF-basic, Carrier-free	12-0005CFR	50 ug, 100 ug, 500 ug, 1 mg
Recombinant Human BCL2	225201	100 ug

Products are for research use only and are not intended for human use. We do not sell to patients.